

OROFACIAL ASSISTIVE DEVICES




Collis Curve Double Action Toothbrush

The toothbrush cleans the outside, inside and chew surfaces at the same time. The bristles of the toothbrush are bent and reach down around the crowns of the teeth. Good for short patience, gap problems and nausea.

Sizes:

- Baby (small brush head)
- Junior (smaller brush head)
- Adult, Medium (large brush head. Medium hard bristles)
- Adult, Soft (large brush head. Soft bristles)
- Adult, Perio (large brush head. Longer brushes)



Proxident Oral Wipe

Proxident oral wipe is extra soft and is used for cleaning the mouth as well as for moisturizing and lubricating dry and irritated mucous membranes. Can be used in conjunction with different fluids and solutions according to the instructions of the dental professional. For single use only. Proxident oral wipe is registered as a medical device and is CE marked.

Sizes:

- Small
- Large



Proxident Bite Support

The bite support is placed between the rows of teeth to prevent the patient from injuring himself on implements/ instruments in the mouth during involuntary biting. Provides relief for the bite while providing support and facilitating for anyone performing oral care or dental treatment. For

frequent use with severe bumps, replacement is recommended after one month. Can be combined with grip ring for a softer and firmer grip.

Material: Polycarbonate plastic, CE marked. Color: white and blue
Cleaning: Disinfection solution or dishwasher max 90° C

Sizes:

- Small
- Large



Bit Spatula Molini mouth opener

Mun-H-Center, in collaboration with Specialized Care Co, Inc., has developed a bit spatula for dental treatment and oral care at home. The bit spatula is placed between the rows of teeth for easier access to the teeth during dental treatment or oral care at home.

The patient cannot injure himself in the event of involuntary biting.

Facilitates keeping your mouth open.

Provides support and facilitation for dental practitioner/oral care practitioner.



Grip enlarger in knurled silicone for better grip

Material: Silicone Cleaning: Machine wash 90 ° C

Sizes:

GF 060 - suitable for the toothpick Easy Pick, TePe

GF 100 - Suitable for double tooth brush, Collis Curve adult, space brush TePe, Glidfloss-Pick dental floss.



Mouth angle expander

The mouth angle holder is used to keep lips and cheeks away during dental treatment, oral care and photography. Provides excellent transparency in treatment and oral care. Available in two sizes, large and small.

Material: Plastic, CE marked Cleaning: Disk disinfectant. Withstands autoclaving at

a maximum of 121 degrees C

Sizes:

Children

Adult

Order at Mun-H-Center

Send an e-mail to mun-h-center@vregion.se including the following information:

Your name/business name

Corporate identity number

Reference

Contact information (e-mail and phone number)

Invoice address

Delivery address

VAT number (Value Added Tax)

Your order (Product name, size and number)

For further help or questions, please contact us at mun-h-center@vregion.se or +46 10 - 441 79 80.